

Nonradioactive lodide Assay kit

Assay kit D05076

- Easy-to-use assay kit
- Nonradioactive waste
- Less than 1 hour of bench time



lodide is essential in the metabolism as it allows temperature regulation, intellectual development for children, muscular development, normal heart function and skeleton growth. lodide transport is the basis for an emerging approach to selective cancer cell destruction.

lodide uptake from blood into thyroid follicular cells is the first step in the biosynthesis of thyroid hormones T3 and T4, known to regulate many essential biological processes.

Thyroid hormones are indispensable for body development. This transport is mediated by NIS (sodium iodide symporter SLC5A5), which is an intrinsic membrane glycoprotein located in the basolateral membrane of thyrocytes. Since the NIS discovery, thorough biochemical analysis has elucidated the mechanism of basolateral lodide transport and revealed the key role of NIS in thyroid diseases such as thyroid cancer, autoimmune diseases, and congenital hypothyroidism.

If the rate is not in normal proportions, some diseases may develop as an underactive thyroid if the rate is too low, or an overactive thyroid if the rate is too high. Other diseases exist such as chronic thyroiditis of Hashimoto, or cancer of the thyroid gland.

lodide deficiency causes many thyroid metabolism disorders; thus, it is important to control lodide rates to prevent all of these diseases.

Endocrine disruptors may affect lodide uptake, leading to thyroid hormone dysregulations.

www.bertin-bioreagent.com www.bertin-corp.com (US webstore)





Nonradioactive Iodide Assay Kit

FOCUS ON THE ASSAY

The present assay is a nonradioactive method for the measurement of lodide.

This lodide assay is based on the oxide reduction reaction:

Cerium(IV) is reduced by arsenic(III). The reduction of yellow (420 nm) cerium(IV) to colorless cerium(III) by arsenic(III) proceeds very slowly but traces of lodide strongly accelerate this reaction with the rate being directly proportional to lodide concentration.

For a given time, the decay is inversely proportional to the lodide concentration in the well.

This method is simple and nonradioactive, and as such it can be used widely.

cat. I⁻ Ce⁴⁺ + ½ As³⁺ → Ce³⁺ + ½ As⁵⁺ Yellow Colorless λmax = 420 nm

TECHNICAL DATA

- Size: 96 wells
- Shelf life: 1-year maximum after production
- Shipping: Wet ice
- Stability: Store at +4°C
- Limit of detection: 0.07 ng/mL

- Sample volume: 100µL
- Standard curve range:100-700 nmoL/L
- Application(s): iodide quantification
- Media application: It is the user's responsibility to check the compatibility of the assay with the study matrix



ntration) nM

RELATED PRODUCTS

20 Hydroxyecdysone ELISA kit

