



Sampling products



**Wide range of pre-analytical tools
now available for your research**



Sampling products flyer

Sample handling has become an essential part in the modern analyst toolbox. In order to get the most reproducible results, it is critical to take care of the samples from the collection to the first step of the analytical procedure.

At early stage after collection, several events may affect the concentration of the targeted compound. In order to fix the concentration as close as possible to the physiological concentrations, Bertin Pharma has developed a **range of solutions ranking from inhibitor cocktails to sampling tubes**.

You will find below our list of **sampling products**:

- | | |
|---|---|
| ▶ Angiotensin II inhibition cocktail
(for up to 5 mL, 10 mL, 50 mL of blood sample) | D05006 |
| Related product
▷ <i>Angiotensin II SPIE-IA kit</i> | <i>A05880</i> |
| ▶ Sampling tubes with BHT (100 ea) | D31007 |
| Related products
▷ <i>8-Isoprostane EIA Kit, oxysterols, ...</i> | <i>516351</i> |
| ▶ Sampling tubes with Indomethacin (100 ea) | D31008 |
| Related products
▷ <i>Thromboxane B2 EIA Kit</i>
▷ <i>2,3-dinor Thromboxane B2 EIA Kit</i>
▷ <i>11-dehydro Thromboxane B2 EIA Kit</i> | <i>519031</i>
<i>519051</i>
<i>519501</i> |
| ▶ Sampling tubes with PHMB (100 ea) | D31009 |
| Related products
▷ <i>Acylated Ghrelin (human) EIA kit</i>
▷ <i>Acylated Ghrelin (rat) EIA kit</i> | <i>A05106</i>
<i>A05117</i> |
| ▶ Sampling tubes with Indomethacin BHT | On request |

Customized sampling tubes with different inhibitors can be prepared upon request.

Please contact us for feasibility and quote.

Your Partner in Innovative Drug Development