



SAMPLING

PRE-ANALYTICAL TOOLS



- 01 Biomarkers
- 02 RTU products for stabilisation
- 03 Validated with our Assays



SAMPLING

ENVIRONMENT & ENDOCRINE DISRUPTORS

Sample handling has become an essential part in the modern analyst toolbox.

In order to get the most reproducible results, it is critical to take care of the samples from the collection to the first step of the analytical procedure.

At early stage after collection, several events may affect the concentration of the target compound.

In order to fix the concentration as close as possible to the physiological concentrations, Bertin Bioreagent has developed a **range of solutions ranking from inhibitor cocktails to sampling tubes**.

**CUSTOMIZED SAMPLING TUBES WITH DIFFERENT INHIBITORS CAN BE PREPARED UPON REQUEST.
PLEASE CONTACT US FOR FEASIBILITY & QUOTE**

SAMPLING PRODUCTS

- | | |
|---|--|
| > Angiotensin II inhibition cocktail (for up to 5mL, 10mL, 50mL of blood sample)
Related product:
- Angiotensin II SPIE-IA kit | Cat No: D05006
Cat No: A05880 |
| > Sampling tubes with BHT (35 or 100 tubes)
Related products:
- 8-Isoprostane EIA kit,
- Oxysterols, ... | Cat No: D31007
Cat No: 516351 |
| > Sampling tubes with Indomethacin (35 or 100 tubes)
Related products:
- Thromboxane B2 EIA kit 519031
- 2,3-dinor Thromboxane B2 EIA kit
- 11-dehydro Thromboxane B2 EIA kit | Cat No: D31008
Cat No: 519031
Cat No: 519051
Cat No: 519501 |
| > Sampling tubes with PHMB (35 or 100 tubes)
Related products:
- Acylated Ghrelin (human) EIA kit
- Acylated Ghrelin (rat) EIA kit | Cat No: D31009
Cat No: A05106
Cat No: A05117 |
| > Sampling tubes with BHT and Indomethacin (35 or 100 tubes), useful for any lipid peroxidation and COX metabolites studies | Cat No: D31010 |