

Steroid Hormone Assay Kits

Researchers have trusted Cayman's collection of steroid hormone assays for decades to study effects on metabolism, inflammation, immune function, reproduction, and stress responses.



Distributed by



Assay Features:

- Highly sensitive detection in the low pg/ml range
- Validated in several biological matrices
- Broad species reactivity

Adrenal Androgens & Androgen Precursor

Item No.	Product Name	Assay Range	Sensitivity (80% B/B ₀)
582701	Testosterone ELISA Kit	3.9-500 pg/ml	6 pg/ml
582751	11-keto Testosterone ELISA Kit	0.78-100 pg/ml	1.3 pg/ml
501740	Testosterone Glucuronide ELISA Kit	4.9-3,000 pg/ml	25 pg/ml
502600	DHEA-S ELISA Kit	31.3-4,000 pg/ml	96 pg/ml

Estrogens

Item No.	Product Name	Assay Range	Sensitivity (80% B/B ₀)
501890	Estradiol ELISA Kit	0.61-10,000 pg/ml	20 pg/ml
501290	Estrone-3-Glucuronide ELISA Kit	2.7-6,000 pg/ml	23.2 pg/ml (overnight), 36.1 pg/ml (rapid)
582261	2-Methoxyestradiol ELISA Kit	13.1-8,000 pg/ml	40 pg/ml

Progestogens

Item No.	Product Name	Assay Range	Sensitivity (80% B/B ₀)
501300	Pregnanediol-3-Glucuronide (PDG) ELISA Kit	0.4-50 ng/ml	0.94 ng/ml
582601	Progesterone ELISA Kit	7.8-1,000 pg/ml	10 pg/ml

Glucocorticoids

Item No.	Product Name	Assay Range	Sensitivity (80% B/B ₀)
501320	Corticosterone ELISA Kit	8.2-5,000 pg/ml	30 pg/ml
500370	Cortisol ELISA Kit (Express)	39.1-5,000 pg/ml	110 pg/ml
501420	11-Oxoetiocholanolone ELISA Kit	0.031-200 ng/ml	0.2 ng/ml

Mineralocorticoids

Item No.	Product Name	Assay Range	Sensitivity (80% B/B ₀)
501090	Aldosterone ELISA Kit	15.6-2,000 pg/ml	30 pg/ml

Vitamin D

Item No.	Product Name	Assay Range	Sensitivity (80% B/B ₀)
501050	Vitamin D ELISA Kit	0.19-25 ng/ml	0.5 ng/ml



Find answers to frequently asked questions about Cayman's steroid hormone kits at www.caymanchem.com/steroidfaq

Explore all of our hormone and receptor assay kits at www.caymanchem.com



Distributed by:
Bertin Bioreagent
Tel: +33 (0)139 306 036
Email: order-life@bertin.fr
Web: www.bertin-bioreagent.com