

Product Information



CGRP (human,rat) Monoclonal Antibody - Clone 83

Cat No: A03481 - 200 µg

General Data

Shipping:	Dry Ice
Formulation:	lyophilized or liquid
Storage Buffer:	Phosphate Buffer
Host:	Mouse
Antigen:	Rat CGRP α -KLH
Clone:	CGRP83
Application(s):	EIA




Image of CGRP (red) from rat sensory neurons labelled with a Monoclonal antibody (A03481 Bertin Pharma) Kindly provided by Neuron Experts

FP/04/24

For research laboratory use only – Not for human diagnostic use.

Buyers agree to purchase the material subject to Purchasing Terms that can be found on our website. Seek appropriate training to safely handle this product under normal conditions. Use the recommended personal protective equipment to prevent chemical exposures.

Bertin Bioreagent does not make any other warranty or representation whatsoever whether expressed or implied, with respect to these products. In no event will Bertin Bioreagent be liable for incidental, consequential or punitive damages.

Contact Bertin Bioreagent

Parc d'Activités du Pas du Lac
10 bis avenue Ampère
78180 Montigny le Bretonneux - France
Tel.: +33 (0)139 306 036

<https://www.bertin-bioreagent.com/pa206/contact-us>
<https://www.bertin-bioreagent.com>