

Product Information



Endothelin (human) Monoclonal Antibody - clone 18

Cat No: A03560 - 200 µg

General Data

Shipping:	Dry Ice
Formulation:	liquid
Storage Buffer:	Phosphate Buffer
Host:	Mouse
Antigen:	Human Endothelin-BSA
Clone:	18
Application(s):	EIA



FP/02/24

For research laboratory use only – Not for human diagnostic use.

Buyers agree to purchase the material subject to Purchasing Terms that can be found on our website. Seek appropriate training to safely handle this product under normal conditions. Use the recommended personal protective equipment to prevent chemical exposures.

Bertin Bioreagent does not make any other warranty or representation whatsoever whether expressed or implied, with respect to these products. In no event will Bertin Bioreagent be liable for incidental, consequential or punitive damages.

Contact Bertin Bioreagent

Parc d'Activités du Pas du Lac
10 bis avenue Ampère
78180 Montigny le Bretonneux - France
Tel.: +33 (0)139 306 036

<https://www.bertin-bioreagent.com/pa206/contact-us>
<https://www.bertin-bioreagent.com>