

Product Information



Goat Anti-Rabbit IgG: DyLight® 650

Cat No: 16560 - 1 mg

General Data

Shipping: wet ice

Formulation: A liquid

Antigen: Rabbit IgG



Application(s): FC, FACS, protein microarrays

Product Overview

Goat anti-Rabbit IgG conjugated to DyLight® 650 Antibody: Goat polyclonal anti-Rabbit IgG Dye: DyLight® 650 Excitation max. λ : 652 nm Emission max. λ : 672 nm Reactivity: Rabbit IgG Immunogen: Rabbit IgG Uses: Flow Cytometry and cell-based assays

FP/19/24

For research laboratory use only – Not for human diagnostic use.

Buyers agree to purchase the material subject to Purchasing Terms that can be found on our website. Seek appropriate training to safely handle this product under normal conditions. Use the recommended personal protective equipment to prevent chemical exposures.

Bertin Bioreagent does not make any other warranty or representation whatsoever whether expressed or implied, with respect to these products. In no event will Bertin Bioreagent be liable for incidental, consequential or punitive damages.

Contact Bertin Bioreagent

Parc d'Activités du Pas du Lac
10 bis avenue Ampère
78180 Montigny le Bretonneux - France
Tel.: +33 (0)139 306 036

<https://www.bertin-bioreagent.com/pa206/contact-us>
<https://www.bertin-bioreagent.com>