

Product Information



Rabbit Anti-Goat IgG: SureLight[®] P3

Cat No: 16698 - 1 mg

General Data

Shipping: wet ice

Formulation: A lyophilized powder



Application(s): FC, FACS, Protein Microarrays

Product Overview

Rabbit anti-Goat IgG, secondary antibody, conjugated to SureLight[®] P3 Antibody: Rabbit polyclonal anti-Goat IgG Isotype: IgG Dye: SureLight[®] P3 Excitation max. λ : 614 nm Emission max. λ : 662 nm Uses: Flow cytometry, cell-based assays, microarrays, and microplate applications Form: Conjugate supplied as a lyophilized powder.

FP/21/24

For research laboratory use only – Not for human diagnostic use.

Buyers agree to purchase the material subject to Purchasing Terms that can be found on our website. Seek appropriate training to safely handle this product under normal conditions. Use the recommended personal protective equipment to prevent chemical exposures.

Bertin Bioreagent does not make any other warranty or representation whatsoever whether expressed or implied, with respect to these products. In no event will Bertin Bioreagent be liable for incidental, consequential or punitive damages.

Contact Bertin Bioreagent

Parc d'Activités du Pas du Lac
10 bis avenue Ampère
78180 Montigny le Bretonneux - France
Tel.: +33 (0)139 306 036

<https://www.bertin-bioreagent.com/pa206/contact-us>
<https://www.bertin-bioreagent.com>