



## Anti-EP<sub>2</sub> Receptor:SureLight™ R-PE

Cat No: 16684 - 50 µg

### General Data

**Shipping:** wet ice

**Formulation:** A lyophilized powder



**Application(s):** FC, FACS, Protein Microarrays

### Product Overview

Anti-EP<sub>2</sub> (Prostaglandin E<sub>2</sub>) Receptor conjugated to SureLight™ R-PE Antibody: Rabbit anti-EP<sub>2</sub> (Prostaglandin E<sub>2</sub>) Receptor  
Reactivity: Reactive with human, murine, and rat EP<sub>2</sub> receptors (other species not tested); Not reactive with EP<sub>1</sub>, EP<sub>3</sub>, and EP<sub>4</sub> receptors. Uses: Flow cytometry, cell-based assays, microarrays, and microplate applications

FP/28/24

**For research laboratory use only – Not for human diagnostic use.**

Buyers agree to purchase the material subject to Purchasing Terms that can be found on our website. Seek appropriate training to safely handle this product under normal conditions. Use the recommended personal protective equipment to prevent chemical exposures.

Bertin Bioreagent does not make any other warranty or representation whatsoever whether expressed or implied, with respect to these products. In no event will Bertin Bioreagent be liable for incidental, consequential or punitive damages.

#### Contact Bertin Bioreagent

Parc d'Activités du Pas du Lac  
10 bis avenue Ampère  
78180 Montigny le Bretonneux - France  
Tel.: +33 (0)139 306 036

<https://www.bertin-bioreagent.com/pa206/contact-us>  
<https://www.bertin-bioreagent.com>