

# Product Information



## PROTEOSAVE SS 1.5mL Microtubes

Cat No: D30006 - 100 ea

### General Data

<b>Shipping:</b>	room temp.
<b>Stability:</b>	Store at RT; shelf life 36 months maximum after production
<b>Formulation:</b>	pack of 100 tubes



### Product Overview

#### Non Protein Adsorption Products - ProteoSaveSS

PROTEOSAVE SS is special surface treatment products with high hydrophilic polymer. Tube or well surfaces are covered with hydroxyl group in high density and it enables to block the non-specific binding of proteins, peptides and cells to the surface dramatically. Sample loss originated from non-specific binding to tubes/plates is prevented effectively. In addition, the property of organic resistance and heat resistance is considerably improved compared with conventional surface treated products.

Reduce the adsorption of proteins, peptides and cells

Resistance to organic solvents, detergents, heat

FP/29/24

**For research laboratory use only – Not for human diagnostic use.**

Buyers agree to purchase the material subject to Purchasing Terms that can be found on our website. Seek appropriate training to safely handle this product under normal conditions. Use the recommended personal protective equipment to prevent chemical exposures.

Bertin Bioreagent does not make any other warranty or representation whatsoever whether expressed or implied, with respect to these products. In no event will Bertin Bioreagent be liable for incidental, consequential or punitive damages.

#### Contact Bertin Bioreagent

Parc d'Activités du Pas du Lac  
10 bis avenue Ampère  
78180 Montigny le Bretonneux - France  
Tel.: +33 (0)139 306 036

<https://www.bertin-bioreagent.com/pa206/contact-us>  
<https://www.bertin-bioreagent.com>