

Product Information



Leukotriene B4 Polyclonal Antibody

Cat No: C03024 - 500 dtn

General Data

Shipping:	wet ice
Formulation:	This vial contains 500 determinations of antibody for a total volume of 1000 µL
Storage Buffer:	Phosphate Buffer
Host:	Rabbit
Antigen:	LTB4-PABSA
Application(s):	RIA
Specificity:	% CROSS REACTIVITY for B/Bo = 50 % Leukotriene B4 100% Leukotriene D4 2,1% Leukotriene E4 8,1% 12-HETE 34%



FP/20/24

For research laboratory use only – Not for human diagnostic use.

Buyers agree to purchase the material subject to Purchasing Terms that can be found on our website. Seek appropriate training to safely handle this product under normal conditions. Use the recommended personal protective equipment to prevent chemical exposures.

Bertin Bioreagent does not make any other warranty or representation whatsoever whether expressed or implied, with respect to these products. In no event will Bertin Bioreagent be liable for incidental, consequential or punitive damages.

Contact Bertin Bioreagent

Parc d'Activités du Pas du Lac
10 bis avenue Ampère
78180 Montigny le Bretonneux - France
Tel.: +33 (0)139 306 036

<https://www.bertin-bioreagent.com/pa206/contact-us>
<https://www.bertin-bioreagent.com>