

Product Information



6,15-diketo Prostaglandin F1 α Polyclonal Antibody

Cat No: C03052 - 500 dtn

General Data

Shipping:	Dry Ice
Formulation:	This vial contains 500 determinations of antibody in 1/2 glycerated liquid for a total volume of 4 μ L.
Host:	Rabbit
Antigen:	6,15-diketo PGF1 α -BSA
Application(s):	RIA
Specificity:	% CROSS REACTIVITY for B/Bo = 50 % E1 < 0,1 K E1 9 DHK E1 0,1 F1 < 0,1 DHF1 50 K F1 3,6 DHKF1 < 0,1 6 keto E1 0,6 6 keto F1 100 6, 15 diketo F1 8,5 E2 < 0,1 DHK E2 0,3 F2 < 0,1 DHKF2 1,2 Thromboxane B2 : (TXB2) < 0,1



FP/16/24

For research laboratory use only – Not for human diagnostic use.

Buyers agree to purchase the material subject to Purchasing Terms that can be found on our website. Seek appropriate training to safely handle this product under normal conditions. Use the recommended personal protective equipment to prevent chemical exposures.

Bertin Bioreagent does not make any other warranty or representation whatsoever whether expressed or implied, with respect to these products. In no event will Bertin Bioreagent be liable for incidental, consequential or punitive damages.

Contact Bertin Bioreagent

Parc d'Activités du Pas du Lac
10 bis avenue Ampère
78180 Montigny le Bretonneux - France
Tel.: +33 (0)139 306 036

<https://www.bertin-bioreagent.com/pa206/contact-us>
<https://www.bertin-bioreagent.com>