

Product Information



TAF10 α Monoclonal Antibody

Cat No: G01076 - 100 μ L

General Data

Shipping: wet ice

Formulation: Liquid Ascite, does not contain any preservative therefore avoid repeat freezing and thawing cycles

Host: Mouse

Antigen: Full length human TAFII30 as a fusion protein

Clone: 2F4

Isotype: IgG1, κ



Application(s): ELISA

Western Blot (endogenous and recombinant proteins)

Immunocytofluorescence

In immunoprecipitation works efficiently, but cannot be eluted with the peptide corresponding to the epitope

Recommended dilutions: 1/500-1/5000

Specificity: Human, does not work on mouse. Recognises the GGPGAGAAAGGTG epitope from aa. 50 to aa. 62

Scientific Literature

Jacq, X., Brou, C., Lutz, Y., Davidson, I., Chambon, P. and Tora, L. (1994) Cell, 79, 107-117

Wieczorek, E., Brand, M., Jacq, X. and Tora, L. (1998) Nature, 393, 187-191

Metzger, D., Scheer, E., Soldatov, A., and Tora, L. (1999) EMBO J., 18, 4823-4834.

FP/19/24

For research laboratory use only – Not for human diagnostic use.

Buyers agree to purchase the material subject to Purchasing Terms that can be found on our website. Seek appropriate training to safely handle this product under normal conditions. Use the recommended personal protective equipment to prevent chemical exposures.

Bertin Bioreagent does not make any other warranty or representation whatsoever whether expressed or implied, with respect to these products. In no event will Bertin Bioreagent be liable for incidental, consequential or punitive damages.

Contact Bertin Bioreagent

Parc d'Activités du Pas du Lac
10 bis avenue Ampère
78180 Montigny le Bretonneux - France
Tel.: +33 (0)139 306 036
<https://www.bertin-bioreagent.com/pa206/contact-us>
<https://www.bertin-bioreagent.com>