



TLF (TBP like factor) Monoclonal Antibody

Cat No: G01105 - 100 µL

General Data

Shipping: wet ice

Formulation: Liquid Ascite, does not contain any preservative therefore avoid repeat freezing and thawing cycles

Host: Mouse

Antigen: Recombinant TLF protein

Clone: 2A1

Isotype: IgG1, κ


Application(s): ELISA

Western Blot

Immunocytochemistry

Recommended dilutions: 1/500-1/5000

Specificity: mouse and human TLF



Expression of TLF in seminiferous tubules of wild-type mice 40x

Scientific Literature

Martianov I, Brancorsini S, Gansmuller A, Parvinen M, Davidson I, Sassone-Corsi P. Distinct functions of TBP and TLF/TRF2 during spermatogenesis: requirement of TLF for heterochromatic chromocenter formation in haploid round spermatids. *Development*. 2002 Feb;129(4):945-55.

Martianov I, Fimia GM, Dierich A, Parvinen M, Sassone-Corsi P, Davidson I. Late arrest of spermiogenesis and germ cell apoptosis in mice lacking the TBP-like TLF/TRF2 gene. *Mol Cell*. 2001 Mar;7(3):509-15.

FP/17/24

For research laboratory use only – Not for human diagnostic use.

Buyers agree to purchase the material subject to Purchasing Terms that can be found on our website. Seek appropriate training to safely handle this product under normal conditions. Use the recommended personal protective equipment to prevent chemical exposures.

Bertin Bioreagent does not make any other warranty or representation whatsoever whether expressed or implied, with respect to these products. In no event will Bertin Bioreagent be liable for incidental, consequential or punitive damages.

Contact Bertin Bioreagent

Parc d'Activités du Pas du Lac
10 bis avenue Ampère
78180 Montigny le Bretonneux - France
Tel.: +33 (0)139 306 036

<https://www.bertin-bioreagent.com/pa206/contact-us>
<https://www.bertin-bioreagent.com>