

# Product Information



## Progesterone receptor $\alpha$ Monoclonal Antibody

Cat No: G01028 - 100  $\mu$ L

### General Data

**Shipping:** wet ice

**Formulation:** Liquid Ascite, does not contain any preservative therefore avoid repeat freezing and thawing cycles

**Host:** Mouse

**Clone:** 6A

**Isotype:** IgG2a,  $\kappa$



**Application(s):** This antibody has been successfully used in immunocytochemistry, (Formaldehyde fixation)  
Immunohistochemistry (frozen section , micro-wave oven treated paraffin sections)  
Elisa tests  
western blot (tested on recombinant protein)  
recommended dilution: 1/1000

**Specificity:** The peptide used as immunogen for this product corresponds to the C terminal sequence of the human PR (end of E region – Hormone binding domain) : This antibody cross reacts with chicken PR and should also work on rabbit PR since the sequence is 100 % con

FP/26/24

For research laboratory use only – Not for human diagnostic use.

Buyers agree to purchase the material subject to Purchasing Terms that can be found on our website. Seek appropriate training to safely handle this product under normal conditions. Use the recommended personal protective equipment to prevent chemical exposures.

Bertin Bioreagent does not make any other warranty or representation whatsoever whether expressed or implied, with respect to these products. In no event will Bertin Bioreagent be liable for incidental, consequential or punitive damages.

### Contact Bertin Bioreagent

Parc d'Activités du Pas du Lac  
10 bis avenue Ampère  
78180 Montigny le Bretonneux - France  
Tel.: +33 (0)139 306 036

<https://www.bertin-bioreagent.com/pa206/contact-us>  
<https://www.bertin-bioreagent.com>