

# Product Information



## Goat Anti-Mouse IgG:DyLight® 650

Cat No: 16559 - 1 mg

### General Data

**Shipping:** wet ice

**Formulation:** A liquid

**Antigen:** Mouse IgG

**Application(s):** FC, FACS, protein microarrays

### Product Overview

Goat anti-Mouse IgG conjugated to DyLight® 650 Antibody: Goat polyclonal anti-Mouse IgG Dye: DyLight® 650 Excitation max.  $\lambda$ : 652 nm Emission max.  $\lambda$ : 672 nm Reactivity: Mouse IgG Immunogen: Mouse IgG Uses: Flow Cytometry and cell-based assays



FP/29/24

For research laboratory use only – Not for human diagnostic use.

Buyers agree to purchase the material subject to Purchasing Terms that can be found on our website. Seek appropriate training to safely handle this product under normal conditions. Use the recommended personal protective equipment to prevent chemical exposures.

Bertin Bioreagent does not make any other warranty or representation whatsoever whether expressed or implied, with respect to these products. In no event will Bertin Bioreagent be liable for incidental, consequential or punitive damages.

#### Contact Bertin Bioreagent

Parc d'Activités du Pas du Lac  
10 bis avenue Ampère  
78180 Montigny le Bretonneux - France  
Tel.: +33 (0)139 306 036

<https://www.bertin-bioreagent.com/pa206/contact-us>  
<https://www.bertin-bioreagent.com>