

Product Information



Goat Anti-Human IgG:SureLight™ APC

Cat No: 17545 - 100 µg

General Data

Shipping: wet ice

Formulation: A lyophilized powder



Application(s): High throughput screening (HTS) and fluorescence resonance energy transfer (FRET) assays, flow cytometry, and microscopy

Product Overview

Goat Anti-Human IgG conjugated to SureLight™ APC Antibody: Goat polyclonal Anti-Human IgG Uses: Optimized for High throughput screening (HTS) and fluorescence resonance energy transfer (FRET) assays, flow cytometry, and microscopy

FP/05/24

For research laboratory use only – Not for human diagnostic use.

Buyers agree to purchase the material subject to Purchasing Terms that can be found on our website. Seek appropriate training to safely handle this product under normal conditions. Use the recommended personal protective equipment to prevent chemical exposures.

Bertin Bioreagent does not make any other warranty or representation whatsoever whether expressed or implied, with respect to these products. In no event will Bertin Bioreagent be liable for incidental, consequential or punitive damages.

Contact Bertin Bioreagent

Parc d'Activités du Pas du Lac
10 bis avenue Ampère
78180 Montigny le Bretonneux - France
Tel.: +33 (0)139 306 036

<https://www.bertin-bioreagent.com/pa206/contact-us>
<https://www.bertin-bioreagent.com>