


Key Benefits

- * Recognized industry leader in Nuclear Receptor research
- * World's largest portfolio of NR Assay Kits (Human & Ortholog)
- * Kits are all-inclusive, cell-based, robust, & ready to use
- * Results are fast, accurate, & easy to analyze

Ready-to-Use NR Assay Systems

- * **Easy Access** - Receive via overnight shipment; store frozen (-80C)
- * **Proven Results** - Robust data, featuring excellent Z' scores and high signal-to-noise ratios
- * **Simple Process** - Go from freezer to data in 24 hours



Each assay kit includes nuclear receptor reporter cells, an assay plate, optimized cell culture media, positive control agonist, and Luciferase Detection Reagent.

- * 3 convenient formats: 3x32-, 96-, & 384-well assay plates
- * Once thawed, cells are ready for immediate use
- * Typically present >95% viability

For all your needs in Nuclear Receptor research, there's only one name you need to know.

distributed by: